

срм-4010 Laboratory Seed Brushing Machine



CLEANING

WORLD LEADER IN SEED PROCESSING EQUIPMENT AND SOLUTIONS



Laboratory Seed Brushing Machine Ø 300 mm CPM-4010

This Laboratory Seed Brushing Machine is used for polishing or threshing (for example seed siliques or umbels) of vegetable, hemp and field (ROW) crop seeds in small batches. It is very suitable for use in breeding stations and commercial processing plants.

The main benefits for you:

- Cleaning of small batches of seeds.
- Threshing of small batches of plant material.
- Easy to operate.
- Compact and easy to clean machine.
- Suitable for a continuous process.

The machine is very suitable for:

Laboratory Seed Brushing Machine Ø 300 mm CPM-4010

Capacity

± 100 grams - 40 kg per hour*

* Indication, depending on crop and settings.

Dimensions

1080 x 700 x 1375 mm (LxWxH)

Noise emission

± 73 dB(A) (measured without product)

Electrical connection

230 V (± 5 %), 1 phase + 0 and ground, 50 / 60 Hz. (Other values are available on request) Installed capacity: 1.6 kW Connection value: 16 A (slow blow)

CE

Including

- The brushes are variable in speed
- Two collecting bins Dust extractor

Part of the CPM Family

уре	Capacity
PM-4000	± 200 - 340 ml per batch
PM-4010	± 100 grams - 40 kg per hour
PM-4020	± 50 - 500 kg per hour
PM-4030	± 80 - 900 kg per hour
PM-4040	± 150 - 1500 kg per hour



Materials

Steel and synthetic materials

Finishing

Powder coated, RAL 9010 (white)

Weight

± 125 kg

Options

- Mantles
- Brushes
- Other inlet and feeding possibilities
- Longer outlets

(One mantle and a set of brushes are necessary for the machine to operate)

Our product groups



EXTRACTION





THRESHING



DISINFECTION

Λ





TREATING





PRIMING

PACKING, COUNTING TRANSPORT & WEIGHING

PORT D



QUALITY CONTROL



HEAD OFFICE Enkhuizen The Netherlands

REGIONAL OFFICE

Salinas The United States of America



Beijing China

Call or mail Seed Processing Holland +31 (0)228 353 000 sales@seedprocessing.nl www.seedprocessing.com

